

## One-Lumen CVC

### CVC KIT: 1-L 16 GA X 50 CM

SKU / Article #: ASK-01651-TRC

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 5**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

One-Lumen CVC SKU / Article #: ASK-01651-TRC	
Quantity	Component Description
1	One Lumen Catheter: 16 Ga. (1.80 mm OD) x 50 cm, Side Holes, Blue FlexTip®
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Syringe: 10 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup

One-Lumen CVC  
 SKU / Article #: ASK-01651-TRC

Quantity	Component Description
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration, with adhesive
1	Needle Holder
1	Scalpel: #11
1	Stopcock: 3-way High-Flow
2	Prep Sponge Swab
6	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 2-0 Silk with Curved Needle
1	Tray: Prep

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.