

## Radial Artery Catheterization

RA CATH KIT: 20 GA X 1-3/4IN (4.45 CM)

SKU / Article #: ASK-04020-UNC

### Highlights

Sharps Safety

### Considerations

Not made with natural rubber latex.  
Includes Lidocaine

Sales Quantity/Case: 10

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Radial Artery Catheterization  
SKU / Article #: ASK-04020-UNC

Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: SafetyGlide™ <sup>1</sup> 18 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep® <sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
1	Wing Clip: Suture
1	SharpsAway® II Locking Disposal Cup
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Needle Holder
1	Safety Scalpel: #11
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>3</sup> 7 cm x 6 cm
1	Suture: 3-0 Silk with Curved Needle
1	Liquid Adhesive: Mastisol® <sup>4</sup>
1	Tubing: 12" (30 cm) Pressure Monitoring





<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup> A registered trademark of 3M Company.

<sup>4</sup> A registered trademark of Ferndale Laboratories, Inc.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.