

## QuickFlash® Radial Artery Catheterization Kit

RA CATH KIT: 20 GA X 1 1/2IN

SKU / Article #: ASK-04220-HHC

### Considerations





Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Representative photo(s): actual product specifications may vary from that shown in the photo

QuickFlash® Radial Artery Catheterization Kit SKU / Article #: ASK-04220-HHC	
Quantity	Component Description
1	QuickFlash® Radial Artery Catheter: 20 Ga. x 1-1/2" (3.81 cm) Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Wing Clip: Suture
1	SharpsAway® II Locking Disposal Cup
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 6.5 cm x 7 cm

## Documents



Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.