

Four-Lumen CVC

CVC KIT: 4-LUMEN 8.5FR X 30CM

SKU / Article #: ASK-14854-KR

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo



Four-Lumen CVC SKU / Article #: ASK-14854-KR	
Quantity	Component Description
1	Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 30 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Slip
1	Syringe: 10 mL Luer-Slip
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm

Four-Lumen CVC
SKU / Article #: ASK-14854-KR

Quantity	Component Description
4	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
2	Medicine Cup: 60 mL
1	SharpsAway® Disposal Cup
1	Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
1	Drape: 30" x 44" (76 cm x 112 cm)
1	Needle Holder
1	Scalpel: #11
4	Stopcock: 4-way
2	Prep Sponge Swab
6	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: OpSite® ¹ 4-3/4" x 4" (12 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle
2	Tray: Prep

¹ A registered trademark of T.J. Smith & Nephew, Limited.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.

© 2025 Teleflex Incorporated. All Rights Reserved.