

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® Two-Lumen CVC Kit

CVC KIT: 2LUMEN 7 FR X 16 CM

SKU / Article #: ASK-22702-CA

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5




Arrowgard Blue[®] Two-Lumen CVC Kit

SKU / Article #: ASK-22702-CA

Quantity	Component Description
1	Arrowgard Blue [®] Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) RW
1	Injection Needle: Eclipse ^{™1} 25 Ga. x 1" (2.54 cm)
1	Pressure Transduction Probe
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow [®] Raulerson Spring-Wire Introduction Syringe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW and 5 mL Luer-Slip Syringe
1	Tissue Dilator
1	SecondSite [™] Adjustable Hub: Fastener
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SharpsAway [®] II Locking Disposal Cup
1	Needle Holder
1	Safety Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Silk with Curved Needle

¹ A trademark of Becton, Dickinson and Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.