

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters

PSI KIT: 9 FR

SKU / Article #: ASK-29903-UHN2

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 5


Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters SKU / Article #: ASK-29903-UHN2	
Quantity	Component Description
1	Arrowgard Blue® Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: SafetyGlide™1 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

Arrowg+ard Blue PSI Kit for use with 7.5 - 8 Fr. Catheters
SKU / Article #: ASK-29903-UHN2

Quantity	Component Description
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Towel
1	Needle Holder
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® ² 10 cm x 12 cm
1	Transducer Cover: 14 cm x 147 cm
1	Suture: 3-0 Nylon with Curved Needle
1	Tubing±: 20" (50 cm) Extension
1	HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.