

## Pressure Injectable Arrowg+ard Blue Plus<sup>®</sup> Three-Lumen CVC Kit

PI CVC KIT: 3L 7 FR X 6IN (16 CM) AGB+

SKU / Article #: ASK-42703-PWO

### Highlights

Antimicrobial: Arrowg+ard Blue Plus\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 3**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

### Pressure Injectable Arrowg+ard Blue Plus<sup>®</sup> Three-Lumen CVC Kit

SKU / Article #: ASK-42703-PWO

Quantity	Component Description
1	Arrowg+ard Blue Plus <sup>®1</sup> Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm, Pressure Injectable
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse <sup>™2</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse <sup>™2</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	SecondSite <sup>™</sup> Adjustable Hub: Fastener
1	SecondSite <sup>™</sup> Adjustable Hub: Catheter Clamp

# Pressure Injectable Arrowgard Blue Plus Three-Lumen CVC Kit

SKU / Article #: ASK-42703-PWO

Quantity	Component Description
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Drape: 77" x 120" Central Line Drape
1	Towel
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Flow Rate Information Card with Injection Log
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>3</sup> 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Braided Silk with Straight Needle
3	Valve: Luer-Activated
1	HemoHopper® Fluid Receptacle






<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of 3M Company.

<sup>4</sup> A registered trademark of ICU Medical, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.