

NO IMAGE
AVAILABLE

Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC Kit/PSI Kit for use with 7 - 7.5 Fr. Catheters

PI CVC KIT: 4L 8.5 FR X 16 CM AGB+

SKU / Article #: ASK-42854-PI1

Highlights

Antimicrobial: Arrowg+ard Blue Plus*

Considerations

*Chlorhexidine

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Plus® Four-Lumen CVC Kit/PSI Kit for use with 7 - 7.5 Fr. Catheters


SKU / Article #: ASK-42854-PI1

Quantity	Component Description
1	Arrowg+ard Blue Plus ^{®1} Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm, Pressure Injectable
1	Arrow-Flex [®] Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow [®] GlideWheel [™] Wire Advancer
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow [®] GlideWheel [™] Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow [®] Raulerson Spring-Wire Introduction Syringe
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
1	Obturator: 8 Fr.
1	Cath-Gard [®] Catheter Contamination Shield: 80 cm with Sterile Tape
1	SecondSite [™] Adjustable Hub: Fastener
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Transducer Cover with Gel and Rubber Bands
4	Valve: MicroCLAVE ^{™2} Neutral Displacement Connector

¹ Licensed under US Patent No. 7,329,412.

² A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.