

## Pressure Injectable Arrowgard Blue Plus® Three-Lumen CVC Kit

CVC KIT: 3 LUMEN 7 FR X 20 CM

SKU / Article #: ASK-45703-PI

### Highlights

---

Antimicrobial: Arrowgard Blue Plus\*  
Sharps Safety

### Considerations

---

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

---

NO IMAGE  
AVAILABLE

*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Pressure Injectable Arrowgard Blue Plus<sup>®1</sup> Three-Lumen CVC Kit



SKU / Article #: ASK-45703-PI

Quantity	Component Description
1	Arrowgard Blue Plus <sup>®1</sup> Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm, Pressure Injectable
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse <sup>TM2</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse <sup>TM2</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	SecondSite <sup>TM</sup> Adjustable Hub: Fastener
1	SecondSite <sup>TM</sup> Adjustable Hub: Catheter Clamp
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	SharpsAway <sup>®</sup> Disposal Cup
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Safety Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.