

## Arrowgard Blue<sup>®</sup> Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

HEMODIALYSIS KIT: 2L 12 FR X 25 CM

SKU / Article #: ASK-46122-VG

### Highlights

Antimicrobial: Arrowgard Blue\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arrowgard Blue <sup>®</sup> Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions	
SKU / Article #: ASK-46122-VG	
Quantity	Component Description
1	Two-Lumen Catheter: 12 Fr. (4.06 mm OD) x 25 cm Radiopaque Polyurethane with Blue FlexTip <sup>®</sup> , Arrowgard <sup>®</sup> Antimicrobial Surface Treatment, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse <sup>™1</sup> 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse <sup>™1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm

# Arrowgard Blue Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions




SKU / Article #: ASK-46122-VG

Quantity	Component Description
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Needle Holder
1	Forceps
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>2</sup> 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Braided Silk with Straight Needle
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper® Fluid Receptacle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.