

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial/Venous Catheterization Set

ART./VENOUS CATH SET:20GA CATH/GA

NDL

SKU / Article #: AV-04020

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 100

Arterial/Venous Catheterization Set SKU / Article #: AV-04020	
Quantity	Component Description
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Wing Clip: Suture

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

