



Nitinol Marked Spring-Wire Guide with  
Arrow® GlideWheel™ Wire Advancer  
Maximum OD per ISO 11070 is 0.89 mm.  
**SPRING WIRE GUIDE COMPONENT**

SKU / Article #: AW-04200-N

**Considerations**


Not made with natural rubber latex.

**Sales Quantity/Case: 25**

Nitinol Marked Spring-Wire Guide with Arrow® GlideWheel™ Wire Advancer Maximum OD per ISO 11070 is 0.89 mm.  
SKU / Article #: AW-04200-N

Quantity	Component Description
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer

**Documents**



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

