

Spring-Wire Guide

SPRING WIRE GUIDE COMPONENT

SKU / Article #: AW-05003

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25

Spring-Wire Guide SKU / Article #: AW-05003	
Quantity	Component Description
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

