

# Marked Spring-Wire Guide with Arrow Advancer

## SPRING WIRE GUIDE COMPONENT

SKU / Article #: AW-14725

### Highlights

CE Mark

### Considerations


Not made with natural rubber latex.

Sales Quantity/Case: 25

Representative photo(s): actual product specifications may vary from that shown in the photo

Marked Spring-Wire Guide with Arrow Advancer SKU / Article #: AW-14725	
Quantity	Component Description
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

