

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC

CVC SET: 3-LUMEN 8.5FR X 16CM

SKU / Article #: BR-12853

Highlights

CE Mark




Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Three-Lumen CVC	
SKU / Article #: BR-12853	
Quantity	Component Description
1	Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11
3	Stopcock: 4-way High-Flow

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.