

## Arrowg+ard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 8.5FR X 16CM

SKU / Article #: BR-22853

### Highlights

---

Antimicrobial: Arrowg+ard Blue\*  
CE Mark

### Considerations

---

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 5**

---




*Representative photo(s): actual product specifications may vary from that shown in the photo*

# Arrowgard Blue® Three-Lumen CVC

SKU / Article #: BR-22853

Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
3	Clamp: Extension Line
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11
3	Stopcock: 4-way High-Flow

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.