

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® Three-Lumen CVC

CVC SET: 3-LUMEN 7FR X 20CM

SKU / Article #: BR-25703

Highlights

Antimicrobial: Arrowgard Blue*
CE Mark




Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrowgard Blue® Three-Lumen CVC SKU / Article #: BR-25703	
Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
3	Dust Cap: Non-Vented
3	Clamp: Extension Line
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Scalpel: #11
3	Stopcock: 4-way High-Flow

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.