



## Radial Artery Catheterization

RA CATH KIT: 20 GA X 1-3/4 IN

SKU / Article #: CA-04020

### Highlights

---

Sharps Safety

### Considerations

---

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

---

*Representative photo(s): actual product specifications may vary from that shown in the photo*


# Radial Artery Catheterization

SKU / Article #: CA-04020

| Quantity | Component Description  |
|----------|--|
| 1        | Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide |
| 1        | Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)  |
| 1        | Syringe: 3 mL Luer-Lock  |
| 1        | SharpsAway® Disposal Cup   |
| 1        | SharpsAway® II Locking Disposal Cup  |
| 1        | Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration   |
| 1        | Towel: 13" x 18" (33 cm x 46 cm)   |
| 1        | Needle Holder  |
| 1        | Safety Scalpel: #11  |
| 1        | Medication Label: 1% Lidocaine   |
| 4        | Gauze Pad: 2" x 2" (5 cm x 5 cm)   |
| 4        | Gauze Pad: 4" x 4" (10 cm x 10 cm)   |
| 1        | Dressing: I.V. Advanced 10 cm x 12 cm  |
| 1        | Surgical Apparel: Bouffant Cap   |
| 1        | Surgical Apparel: Mask with Eye Shield   |
| 1        | Suture: 3-0 Silk with Curved Needle  |
| 1        | Tubing±: 6" Extension with One-Way Stopcock  |

<sup>1</sup> A trademark of Becton, Dickinson and Company.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.