



Central Venous Catheterization Set

CVC SET: 14 GA X 6"

SKU / Article #: CA-04706

Contraindications & Considerations


This product contains dry natural rubber.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Central Venous Catheterization Set SKU / Article #: CA-04706	
Quantity	Component Description
1	One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 16 cm, Side Holes
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	T-Port Connector
1	Scalpel: #11

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.