



## Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

PSI KIT: 8.5 FR

SKU / Article #: CA-09886

### Highlights

CE Mark

### Considerations

Caution: This product contains dry natural rubber.

Sales Quantity/Case: 10



*Representative photo(s): actual product specifications may vary from that shown in the photo*

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

SKU / Article #: CA-09886

Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
1	Pressure Transduction Probe
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Scalpel: #11
1	Stopcock: 3-way High-Flow
2	Prep Sponge Swab
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.