

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue® Two-Lumen CVC

CVC KIT: 2-LUMEN 7 FR X 20 CM

SKU / Article #: CA-27702

Highlights

Antimicrobial: Arrowg+ard Blue*

Considerations


*Chlorhexidine

Not made with natural rubber latex.

Sales Quantity/Case: 5

Arrowg+ard Blue® Two-Lumen CVC SKU / Article #: CA-27702	
Quantity	Component Description
1	Arrowg+ard Blue® Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® Disposal Cup
1	Drape: 21" x 36" (53 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.