

NO IMAGE  
AVAILABLE

## Arterial Catheterization Kit

### ARTERIAL CATH KIT: 20 GA

SKU / Article #: CAN-00820-EDM

#### Highlights

Sharps Safety

#### Considerations

Not made with natural rubber latex.


**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

<b>Arterial Catheterization Kit</b> <b>SKU / Article #: CAN-00820-EDM</b>	
Quantity	Component Description
1	Catheter: 20 Ga. x 8 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: Needle-Pro <sup>®1</sup> 25 Ga. x 5/8" (1.60 cm) UTW
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SharpsAway <sup>®</sup> Disposal Cup
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	O.R. Towels (4 Per Pack)
1	Needle Holder
1	Scissors
6	Gauze Pad: 2" x 2" (5 cm x 5 cm)
1	Dressing: Tegaderm <sup>®2</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover: 14 cm x 147 cm
1	Suture: 3-0 Nylon with Curved Needle

<sup>1</sup> A registered trademark of Portex, Inc.  
<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

