

NO IMAGE
AVAILABLE

Access Tray with TipTracker™ Stylet for use with 5 Fr. PICC

PICC ACCESS TRAY

SKU / Article #: CAN-04001-AX

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 3

Representative photo(s): actual product specifications may vary from that shown in the photo

Access Tray with TipTracker™ Stylet for use with 5 Fr. PICC

SKU / Article #: CAN-04001-AX


Quantity	Component Description
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	VPS TipTracker™ Stylet with T-Port Connector: 0.015" (0.38 mm) dia. x 29" (74 cm)
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer

Access Tray with TipTracker Stylet for use with 5 Fr. PICC
SKU / Article #: CAN-04001-AX

Quantity	Component Description
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
2	Towel
1	Needle Holder
1	Sterile Procedure Sign
2	Paper Tape Measure
2	Tourniquet
1	Foam Electrodes (3 per pack)
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 10 cm x 15.5 cm
1	Transducer Cover with Gel and Rubber Bands
1	Remote Cover
1	Surgical Tape: Transpore™ ² 1" x 54" Roll
1	Valve: MaxZero™ ³ Needleless Connector
1	T-Piece Cover

¹ A trademark of Becton, Dickinson and Company.
² A trademark of 3M Company.
³ A trademark of CareFusion or one of its subsidiaries.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.