

NO IMAGE
AVAILABLE

Arterial Catheterization Kit

RA CATH KIT: 20 GA X 4.45 CM

SKU / Article #: CAN-04020-NSH

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo


Arterial Catheterization Kit
SKU / Article #: CAN-04020-NSH

Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
1	Towel: 13" x 18" (33 cm x 46 cm)
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 10 cm x 12 cm
1	Transducer Cover: 14 cm x 147 cm
1	Dressing: STATLOCK®2 Arterial Select Catheter Stabilization Device 8.9 cm x 3.8 cm
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of C. R. Bard, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for

Use available under the Documents section on this page, where an IFU is available.