

NO IMAGE
AVAILABLE

PSI Kit for use with 7 - 7.5 Fr. Catheters

PSI KIT: 8.5 FR

SKU / Article #: CAN-09803

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo


PSI Kit for use with 7 - 7.5 Fr. Catheters

SKU / Article #: CAN-09803

Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: Eclipse™ ¹ 23 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
1	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Drape: 30" x 36" (76 cm x 91 cm) with 3" (7.6 cm) fenestration, with adhesive
1	Safety Scalpel: #11
1	Stopcock: 3-way High-Flow
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Transducer Cover: 14 cm x 147 cm
1	Suture: 3-0 Nylon with Straight Needle
1	Tubing±: 20" (50 cm) Extension

¹ A trademark of Becton, Dickinson and Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.