

## Arrowgard Blue® MAC™ Two-Lumen Central Venous Access Kit for use with 7.5 - 8 Fr. Catheters

MAC KIT: 2L 9 FR DISTAL X 11.5 CM AGB

SKU / Article #: CAN-21242-XCN1A

### Highlights

Antimicrobial: Arrowgard Blue\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arrowgard Blue® MAC™ Two-Lumen Central Venous Access Kit for use with 7.5 - 8 Fr. Catheters SKU / Article #: CAN-21242-XCN1A	
Quantity	Component Description
1	Two-Lumen Central Venous Access Device: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Arrowgard® Antimicrobial Surface Treatment <sup>1</sup> , Integral Hemostasis Valve, Extension Line Clamps
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse™ <sup>2</sup> 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse™ <sup>2</sup> 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Tissue Dilator
1	Obturator: 8 Fr.

Arrowgard Blue MAC Two-Lumen Central Venous Access Kit for use with 7.5 - 8 Fr. Catheters  
SKU / Article #: CAN-21242-XCN1A




Quantity	Component Description
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	Dust Cap: Vented
1	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: Tegaderm® <sup>3</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover: 14 cm x 147 cm
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Tubing±: 20" (50 cm) Extension
1	HemoHopper® Fluid Receptacle

<sup>1</sup> Licensed under US Patent No. 6,706,024.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of 3M Company.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.