

Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters

PSI KIT: 8.5 FR X 10 CM ANTIMICROBIAL

SKU / Article #: CAN-29803-XCN1A

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
±Precaution: Contains Phthalate: DEHP

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowgard Blue® PSI Kit for use with 7.5 - 8 Fr. Catheters	
SKU / Article #: CAN-29803-XCN1A	
Quantity	Component Description
1	Arrowgard Blue® Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse™1 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	Dust Cap: Non-Vented




Arrowgard Blue PSI Kit for use with 7.5 - 8 Fr. Catheters
SKU / Article #: CAN-29803-XCN1A

Quantity	Component Description
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: Tegaderm® ² CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover: 14 cm x 147 cm
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Tubing±: 20" (50 cm) Extension
1	HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.