

## Pressure Injectable Arrowg+ard Blue Plus<sup>®</sup> Four-Lumen CVC Kit

PI CVC KIT: 4L 8.5 FR X 16 CM AGB+

SKU / Article #: CAN-42854-XPCN1A

### Highlights

Antimicrobial: Arrowg+ard Blue Plus\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
±Precaution: Contains Phthalate: DEHP

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

| Pressure Injectable Arrowg+ard Blue Plus <sup>®</sup> Four-Lumen CVC Kit |  |
|--|--|
| SKU / Article #: CAN-42854-XPCN1A  |  |
| Quantity   | Component Description  |
| 1  | Arrowg+ard Blue Plus <sup>®1</sup> Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 16 cm, Pressure Injectable   |
| 1  | Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow <sup>®</sup> GlideWheel <sup>™</sup> Wire Advancer |
| 1  | Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle  |
| 1  | Injection Needle: Eclipse <sup>™2</sup> 25 Ga. x 1-1/2" (3.81 cm)  |
| 1  | Injection Needle: Eclipse <sup>™2</sup> 23 Ga. x 1-1/2" (3.81 cm)  |
| 1  | Pressure Transduction Probe  |
| 1  | Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow <sup>®</sup> Raulerson Spring-Wire Introduction Syringe  |
| 1  | Syringe: 3 mL Luer-Lock  |
| 1  | Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm  |
| 1  | SecondSite <sup>™</sup> Adjustable Hub: Fastener   |

Pressure Injectable Arrowgard Blue Plus Four-Lumen CVC Kit

SKU / Article #: CAN-42854-XPCN1A




| Quantity | Component Description   |
|----------|---|
| 1        | SecondSite™ Adjustable Hub: Catheter Clamp                              |
| 1        | SharpsAway® II Locking Disposal Cup                                     |
| 1        | SharpsAway® Disposal Cup  |
| 1        | Maximal Barrier Drape™ with 4" fenestration                             |
| 1        | Towel   |
| 1        | Needle Holder   |
| 1        | Forceps   |
| 1        | Safety Scalpel: #11   |
| 1        | Checklist/CLIP Sheet  |
| 1        | Sterile Procedure Sign  |
| 2        | Gauze Pad: 2" x 2" (5 cm x 5 cm)  |
| 10       | Gauze Pad: 4" x 4" (10 cm x 10 cm)                                      |
| 1        | Surgical Apparel: Impervious Gown                                       |
| 1        | Dressing: Tegaderm® <sup>3</sup> CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm) |
| 1        | Transducer Cover: 14 cm x 147 cm  |
| 1        | Surgical Apparel: Bouffant Cap  |
| 1        | Surgical Apparel: Mask with Eye Shield                                  |
| 1        | Suture: 3-0 Silk with Curved Needle                                     |
| 1        | Tubing±: 20" (50 cm) Extension  |
| 1        | HemoHopper® Fluid Receptacle  |

<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of 3M Company.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.