

NO IMAGE
AVAILABLE

Arrowgard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions

LBCVC AGB KIT: 3-L 12 FR X 20 CM

SKU / Article #: CDA-25123-XCN1A

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo


Arrowgard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions SKU / Article #: CDA-25123-XCN1A	
Quantity	Component Description
1	Arrowgard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: Eclipse™1 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Injection Needle: Eclipse™1 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 10 mL Luer-Lock
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
1	SharpsAway® II Locking Disposal Cup

Arrowg+ard Blue Large-Bore Three-Lumen CVC for High Volume Infusions SKU / Article #: CDA-25123-XCN1A	
Quantity	Component Description
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Needle Holder
1	Forceps
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: Tegaderm® ² CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)
1	Transducer Cover: 14 cm x 147 cm
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.