

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC

OUS PI PICC KIT: 1-L 4.5FR X 40CM CG+

SKU / Article #: CDA-44041-HPK1A

Highlights

- Antimicrobial: Arrowg+ard Blue Advance*
- Antithrombogenic: Arrowg+ard Blue Advance*
- Sharps Safety
- CE Mark

Considerations

- *Chlorhexidine
- Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

| Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC | |
|---|--|
| SKU / Article #: CDA-44041-HPK1A | |
| Quantity | Component Description |
| 1 | One Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 4.5 Fr. (1.58 mm OD) x 40 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire |
| 1 | GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator |
| 1 | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer |
| 1 | Injection Needle: Eclipse™ ² 25 Ga. x 1" (2.54 cm) |
| 1 | Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW |
| 1 | Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW |
| 1 | Syringe: 10 mL Luer-Lock |
| 1 | Syringe: 3 mL Luer-Lock |
| 1 | Dust Cap: Non-Vented |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC

SKU / Article #: CDA-44041-HPK1A

| Quantity | Component Description |
|----------|--|
| 1 | SecondSite™ Adjustable Hub: Fastener |
| 1 | SharpsAway® II Locking Disposal Cup |
| 1 | SharpsAway® Disposal Cup |
| 1 | Catheter Trimmer |
| 1 | Towel: 24" x 36" (61 cm x 91 cm) |
| 1 | Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations |
| 1 | Towel |
| 1 | Safety Scalpel: #11 |
| 1 | Patient ID Card |
| 1 | Patient Information Booklet |
| 2 | Paper Tape Measure |
| 1 | Tourniquet |
| 10 | Gauze Pad: 2" x 2" (5 cm x 5 cm) |
| 10 | Gauze Pad: 4" x 4" (10 cm x 10 cm) |
| 1 | Surgical Apparel: Impervious Gown |
| 1 | Dressing: I.V. Advanced 8.5 cm x 11.5 cm |
| 1 | Surgical Apparel: Mask with Eye Shield |
| 1 | Dressing: STATLOCK® ³ PICC Plus Catheter Stabilization Device |
| 1 | Surgical Apparel: Bouffant Cap |
| 1 | Surgical Apparel: Mask |
| 1 | Tape: Steri-Strip® ⁴ |
| 1 | HemoHopper® Fluid Receptacle |





¹ Licensed under US Patent No. 7,329,412.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of C. R. Bard, Inc.

⁴ A registered trademark of 3M Company.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.