

## Pressure Injectable Arrowgard Blue Advance® Three-Lumen PICC

PI CG+ PICC KIT: 3-L 6 FR X 50 CM

SKU / Article #: CDA-45063-HPK1A

### Highlights

- Antimicrobial: Arrowgard Blue Advance\*
- Antithrombogenic: Arrowgard Blue Advance\*
- Sharps Safety
- CE Mark

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Pressure Injectable Arrowgard Blue Advance® Three-Lumen PICC	
SKU / Article #: CDA-45063-HPK1A	
Quantity	Component Description
1	Three Lumen TaperFree® Catheter with Arrowgard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 6 Fr. (2.11 mm OD) x 50 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire
1	GlideThru™ Peel-Away Sheath: 6 Fr. x 2-3/4" (7 cm) Radiopaque over 6 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: Eclipse™ <sup>2</sup> 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
3	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener

Pressure Injectable Arrowgard Blue Advance Three-Lumen PICC  
 SKU / Article #: CDA-45063-HPK1A

Quantity	Component Description
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>3</sup> 15.5 cm x 10 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: STATLOCK® <sup>4</sup> Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>3</sup>
1	HemoHopper® Fluid Receptacle





<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of 3M Company.

<sup>4</sup> A registered trademark of C. R. Bard, Inc.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.