

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

AGBA PICC/DELTA KIT: 1-L 4.5 FR X 55 CM

SKU / Article #: CDA-45541-VPS2

Highlights

- Antimicrobial: Arrowg+ard Blue Advance*
- Antithrombogenic: Arrowg+ard Blue Advance*
- Sharps Safety
- Pre-loaded with VPS Stylet

Considerations

- *Chlorhexidine
- Not made with natural rubber latex.

Sales Quantity/Case: 3

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

SKU / Article #: CDA-45541-VPS2

Quantity	Component Description
1	One Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 4.5 Fr. (1.58 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®, Arrow® VPS Precision™ Stylet, Tuohy-Borst Adapter, Marking Device, and Console Pointers
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ ² 25 Ga. x 5/8" (1.60 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet

SKU / Article #: CDA-45541-VPS2

Quantity	Component Description
1	Syringe: 3 mL Luer-Lock
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® ³ PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Surgical Tape: Transpore™ ⁴ 1" x 54" Roll

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet
 SKU / Article #: CDA-45541-VPS2

Quantity	Component Description
1	HemoHopper® Fluid Receptacle





¹ Licensed under US Patent No. 7,329,412.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of C. R. Bard, Inc.

⁴ A trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.