

## Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC

OUS PI PICC KIT: 2-L 5.5FR X 55CM CG+

SKU / Article #: CDA-45552-HPK1A

### Highlights

- Antimicrobial: Arrowg+ard Blue Advance\*
- Antithrombogenic: Arrowg+ard Blue Advance\*
- Sharps Safety
- CE Mark

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.

**Sales Quantity/Case: 5**

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC	
SKU / Article #: CDA-45552-HPK1A	
Quantity	Component Description
1	Two Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, T-Port Connector, Blue FlexTip® and Placement Wire
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: Eclipse™ <sup>2</sup> 25 Ga. x 1" (2.54 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC

SKU / Article #: CDA-45552-HPK1A

Quantity	Component Description
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Patient ID Card
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>3</sup> 15.5 cm x 10 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: STATLOCK® <sup>4</sup> Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>3</sup>
1	HemoHopper® Fluid Receptacle




<sup>1</sup> Licensed under US Patent No. 7,329,412.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

<sup>3</sup> A registered trademark of 3M Company.

<sup>4</sup> A registered trademark of C. R. Bard, Inc.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.