



Representative photo(s): actual product specifications may vary from that shown in the photo

One-Lumen Midline

MIDLINE CATHETER KIT: 4 FR X 20 CM

SKU / Article #: CDC-02041-MK1A

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 5

One-Lumen Midline SKU / Article #: CDC-02041-MK1A	
Quantity	Component Description
1	Midline Catheter, Trimmable: 4 Fr. (1.4 mm OD) x 20 cm, T-Port Connector, Contamination Guard and Placement Wire
1	GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
1	Injection Needle: SafetyGlide™1 25 Ga. x 5/8" (1.60 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep®2 One-Step Solution with Hi-Lite Orange™ Tint
1	Dust Cap: Non-Vented

One-Lumen Midline
SKU / Article #: CDC-02041-MK1A

Quantity	Component Description
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
5	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH® ³
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® ⁴ PICC Plus Catheter Stabilization Device
1	Tape: Steri-Strip® ⁵
1	HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.





² A registered trademark of CareFusion or one of its subsidiaries.

³ A registered trademark of Johnson & Johnson Corporation.

⁴ A registered trademark of C. R. Bard, Inc.

⁵ A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.