

Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions

LBCVC KIT: 3L 12 FR X 16 CM AGB

SKU / Article #: CDC-22123-1A

Highlights

Antimicrobial: Arrowg+ard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 5

Arrowg+ard Blue® Large-Bore Three-Lumen CVC for High Volume Infusions

SKU / Article #: CDC-22123-1A


Quantity	Component Description
1	Arrowg+ard Blue® Three-Lumen Catheter: 12 Fr. (4.1 mm OD) x 16 cm
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: Eclipse™ 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: SafetyGlide™ 23 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm

Arrowg+ard Blue Large-Bore Three-Lumen CVC for High Volume Infusions
SKU / Article #: CDC-22123-1A

Quantity	Component Description
1	Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep® ² One-Step Solution with Hi-Lite Orange™ Tint
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.