

## Small Vein Arrowgard Blue® Two-Lumen CVC

2-L CVC KIT: 4 FR X 5 CM

SKU / Article #: CDC-22402-X1A

### Highlights

Antimicrobial: Arrowgard Blue\*  
Sharps Safety

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
Includes Lidocaine

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Small Vein Arrowgard Blue® Two-Lumen CVC	
SKU / Article #: CDC-22402-X1A	
Quantity	Component Description
1	Arrowgard Blue® Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 5 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1.16" (30 mm) Insyte™ Autoguard™ Radiopaque over 22 Ga. Introducer Needle
1	Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW and 5 mL Luer-Slip Syringe
1	Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
2	10 mL Luer-Lock Pre-Filled Saline Syringe

Small Vein Arrowgard Blue Two-Lumen CVC

SKU / Article #: CDC-22402-X1A

Quantity	Component Description
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®2</sup> One-Step Solution with Hi-Lite Orange™ Tint
2	Dust Cap: Non-Vented
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	SharpsAway <sup>®</sup> Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Towel
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm <sup>®3</sup> 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: BIOPATCH <sup>®4</sup>
1	Dressing: STATLOCK <sup>®5</sup> CV Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper <sup>®</sup> Fluid Receptacle

<sup>1</sup> A trademark of Becton, Dickinson and Company.


<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup> A registered trademark of 3M Company.



<sup>4</sup> A registered trademark of Johnson & Johnson Corporation.

<sup>5</sup> A registered trademark of C. R. Bard, Inc.

# Documents



Model	Length	Gauge	Quantity
SafetyGlide™	1.5"	25G	100
SafetyGlide™	1.5"	27G	100
SafetyGlide™	1.5"	30G	100
SafetyGlide™	1.5"	32G	100
SafetyGlide™	1.5"	34G	100
SafetyGlide™	1.5"	36G	100
SafetyGlide™	1.5"	38G	100
SafetyGlide™	1.5"	40G	100
SafetyGlide™	1.5"	42G	100
SafetyGlide™	1.5"	44G	100
SafetyGlide™	1.5"	46G	100
SafetyGlide™	1.5"	48G	100
SafetyGlide™	1.5"	50G	100



Model	Length	Gauge	Quantity
SafetyGlide™	1.5"	25G	100
SafetyGlide™	1.5"	27G	100
SafetyGlide™	1.5"	30G	100
SafetyGlide™	1.5"	32G	100
SafetyGlide™	1.5"	34G	100
SafetyGlide™	1.5"	36G	100
SafetyGlide™	1.5"	38G	100
SafetyGlide™	1.5"	40G	100
SafetyGlide™	1.5"	42G	100
SafetyGlide™	1.5"	44G	100
SafetyGlide™	1.5"	46G	100
SafetyGlide™	1.5"	48G	100
SafetyGlide™	1.5"	50G	100

SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.