

Small Vein Arrowgard Blue® Two-Lumen CVC

2-L CVC KIT: 4 FR X13 CM

SKU / Article #: CDC-24402-X1A

Highlights

Antimicrobial: Arrowgard Blue*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Small Vein Arrowgard Blue® Two-Lumen CVC	
SKU / Article #: CDC-24402-X1A	
Quantity	Component Description
1	Arrowgard Blue® Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
1	Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1.16" (30 mm) Insyte™1 Autoguard™ Radiopaque over 22 Ga. Introducer Needle
1	Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW and 5 mL Luer-Slip Syringe
1	Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
1	Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm

Small Vein Arrowgard Blue Two-Lumen CVC
 SKU / Article #: CDC-24402-X1A

Quantity	Component Description
2	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®2} One-Step Solution with Hi-Lite Orange [™] Tint
2	Dust Cap: Non-Vented
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SecondSite [™] Adjustable Hub: Fastener
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Towel
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm ^{®3} 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: BIOPATCH ^{®4}
1	Dressing: STATLOCK ^{®5} CV Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper [®] Fluid Receptacle


¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of CareFusion or one of its subsidiaries.

³ A registered trademark of 3M Company.

⁴ A registered trademark of Johnson & Johnson Corporation.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.