



*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Arrowgard Blue® PSI for use with 7.5 - 8 Fr. Catheters

### PSI KIT: 9 FR X 10 CM ANTIMICROBIAL

SKU / Article #: CDC-29903-XN1A

#### Highlights

Antimicrobial: Arrowgard Blue\*  
Sharps Safety

#### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.  
Includes Lidocaine  
±Precaution: Contains Phthalate: DEHP

#### Sales Quantity/Case: 5

### Arrowgard Blue® PSI for use with 7.5 - 8 Fr. Catheters

SKU / Article #: CDC-29903-XN1A

Quantity	Component Description
1	Arrowgard Blue® Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Nitinol, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow® GlideWheel™ Wire Advancer
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: SafetyGlide™1 25 Ga. x 1" (2.54 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Pressure Transduction Probe
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
1	Injection Needle: SafetyGlide™1 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
1	Syringe: 3 mL Luer-Lock
2	10 mL Luer-Lock Pre-Filled Saline Syringe

Arrowgard Blue PSI for use with 7.5 - 8 Fr. Catheters

SKU / Article #: CDC-29903-XN1A







Quantity	Component Description
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep <sup>®2</sup> One-Step Solution with Hi-Lite Orange™ Tint
1	Obturator: 8 Fr.
1	Cath-Gard <sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	Dust Cap: Non-Vented
1	SharpsAway <sup>®</sup> Disposal Cup
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Needle Holder
1	Safety Scalpel: #11
1	Stopcock: 4-way High-Flow
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Medication Label: 1% Lidocaine
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
5	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH <sup>®3</sup>
1	Transducer Cover: 14 cm x 14.7 cm
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Suture: 3-0 Silk with Curved Needle
1	Tubing±: 20" (50 cm) Extension
1	HemoHopper <sup>®</sup> Fluid Receptacle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup> A registered trademark of Johnson & Johnson Corporation.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2024 Teleflex Incorporated. All Rights Reserved.