

## Pressure Injectable Two-Lumen Midline

### PI MIDLINE 2-L: 5FR X 20CM

SKU / Article #: CDC-32052-MPKB

### Highlights

Sharps Safety


### Considerations

Not made with natural rubber latex.

Includes Lidocaine

###MR Safe

**Sales Quantity/Case: 5**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

Pressure Injectable Two-Lumen Midline SKU / Article #: CDC-32052-MPKB	
Quantity	Component Description
1	Two-Lumen TaperFree® Catheter###: 5 Fr.(1.81 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with Extension Line Clamps and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 5 Fr. x 2-3/4" (7 cm) Radiopaque over 5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™1 25 Ga. x 5/8" (1.60 cm) RW
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Lock
2	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep®2 One-Step Solution
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep®2 One-Step Solution with Hi-Lite Orange™ Tint

## Pressure Injectable Two-Lumen Midline

SKU / Article #: CDC-32052-MPKB

Quantity	Component Description
2	Dust Cap: Non-Vented
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Overbed Sign
1	Medication Label: 1% Lidocaine
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH® <sup>3</sup>
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK® <sup>4</sup> PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Tape: Steri-Strip® <sup>5</sup>

<sup>1</sup> A trademark of Becton, Dickinson and Company.







<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup> A registered trademark of Johnson & Johnson Corporation.

<sup>4</sup> A registered trademark of C. R. Bard, Inc.

<sup>5</sup> A registered trademark of 3M Company.

# Documents



[SCROLL TO VIEW MORE →](#)

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.