

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline

AGBA PI MIDLINE 2-L: 5.5FR X 15CM

SKU / Article #: CDC-41552-MPKB

Highlights

Antimicrobial: Arrowg+ard Blue Advance*

Antithrombogenic: Arrowg+ard Blue Advance*

Sharps Safety

Considerations

*Chlorhexidine

Not made with natural rubber latex.

Includes Lidocaine

†††MR Safe

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen Midline SKU / Article #: CDC-41552-MPKB	
Quantity	Component Description
1	Two-Lumen TaperFree® Catheter††† with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 15 cm Pressure Injectable with Blue FlexTip® Feature and Contamination Guard
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ ² 25 Ga. x 5/8" (1.60 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Syringe: 5 mL Luer-Lock
2	10 mL Luer-Lock Pre-Filled Saline Syringe


Pressure Injectable Arrowg+ard Blue Advance Two-Lumen Midline
SKU / Article #: CDC-41552-MPKB

Quantity	Component Description
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep®3 One-Step Solution
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep®3 One-Step Solution with Hi-Lite Orange™ Tint
2	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with two 13 cm (5") x 15 cm (6") fenestrations
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Overbed Sign
1	Medication Label: 1% Lidocaine
1	Patient Information Booklet
2	Paper Tape Measure
1	Tourniquet
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH®4
1	Transducer Cover with Gel and Rubber Bands
1	Surgical Apparel: Mask with Eye Shield
1	Dressing: STATLOCK®5 PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap

Pressure Injectable Arrowg+ard Blue Advance Two-Lumen Midline SKU / Article #: CDC-41552-MPKB	
Quantity	Component Description
1	Surgical Apparel: Mask
1	Tape: Steri-Strip®6

¹ Licensed under US Patent No. 7,329,412 and 11,219,706.
² A trademark of Becton, Dickinson and Company.
³ A registered trademark of CareFusion or one of its subsidiaries.
⁴ A registered trademark of Johnson & Johnson Corporation.
⁵ A registered trademark of C. R. Bard, Inc.
⁶ A registered trademark of 3M Company.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.