

NO IMAGE
AVAILABLE

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen JACC

PI JACC KIT: 1-L 4.5 FR X 25 CM CG+

SKU / Article #: CDC-42541-JX1A

Highlights

Antimicrobial: Arrowg+ard Blue Advance*
Antithrombogenic: Arrowg+ard Blue Advance*
Sharps Safety

Considerations

*Chlorhexidine
Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 3

Pressure Injectable Arrowg+ard Blue Advance® One-Lumen JACC SKU / Article #: CDC-42541-JX1A	
Quantity	Component Description
1	One-Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 4.5 Fr. (1.55 mm OD) x 25 cm, Pressure Injectable, Blue FlexTip®
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Catheter: 20 Ga. x 1-3/4" (4.50 cm) over 22 Ga. Safety Introducer Needle
1	Injection Needle: Eclipse™ ² 25 Ga. x 1-1/2" (3.81 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 5 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip

Pressure Injectable Arrowguard Blue Advance One-Lumen JACC

SKU / Article #: CDC-42541-JX1A

Quantity	Component Description
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
2	3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint
1	Dust Cap: Non-Vented
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	SecondSite™ Adjustable Hub: Fastener
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape™ with 4" fenestration
1	Towel
1	Safety Scalpel: #11
1	Checklist/CLIP Sheet
1	Sterile Procedure Sign
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
2	Paper Tape Measure
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH®4
1	Dressing: STATLOCK®5 CV Plus Catheter Stabilization Device
1	Surgical Apparel: Mask with Eye Shield
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Transducer Cover: 5" x 48" (13 cm x 122 cm) with Sterile Gel


Pressure Injectable Arrowg+ard Blue Advance One-Lumen JACC

SKU / Article #: CDC-42541-JX1A

Quantity	Component Description
1	HemoHopper® Fluid Receptacle

¹ Licensed under US Patent No. 7,329,412.
² A trademark of Becton, Dickinson and Company.
³ A registered trademark of CareFusion or one of its subsidiaries.
⁴ A registered trademark of Johnson & Johnson Corporation.
⁵ A registered trademark of C. R. Bard, Inc.

Documents



SCROLL TO VIEW MORE ➔

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.