

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

**AGBA PICC/DELTA KIT: 1-L 4.5 FR X 55 CM**

SKU / Article #: CDC-45541-VPS2

### Highlights

- Antimicrobial: Arrowg+ard Blue Advance\*
- Antithrombogenic: Arrowg+ard Blue Advance\*
- Sharps Safety
- Pre-loaded with VPS Stylet

### Considerations

- \*Chlorhexidine
- Not made with natural rubber latex.
- Includes Lidocaine

**Sales Quantity/Case: 3**

### Pressure Injectable Arrowg+ard Blue Advance® One-Lumen PICC pre-loaded with Arrow® VPS Precision™ Stylet

SKU / Article #: CDC-45541-VPS2

Quantity	Component Description
1	One Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection <sup>1</sup> : 4.5 Fr. (1.58 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®, Arrow® VPS Precision™ Stylet, Tuohy-Borst Adapter, Marking Device, and Console Pointers
1	GlideThru™ Peel-Away Sheath: 4.5 Fr. x 2-3/4" (7 cm) Radiopaque over 4.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ <sup>2</sup> 25 Ga. x 5/8" (1.60 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

Pressure Injectable Arrow+ard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet

SKU / Article #: CDC-45541-VPS2


Quantity	Component Description
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA ChloroPrep <sup>®3</sup> One-Step Solution
1	3 mL Applicator 2% CHG and 70% IPA ChloroPrep <sup>®3</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint
1	Dust Cap: Non-Vented
1	SecondSite <sup>™</sup> Adjustable Hub: Catheter Clamp
1	SecondSite <sup>™</sup> Adjustable Hub: Fastener
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	SharpsAway <sup>®</sup> Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape <sup>™</sup> with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)

Pressure Injectable Arrowgard Blue Advance One-Lumen PICC pre-loaded with Arrow VPS Precision Stylet  
 SKU / Article #: CDC-45541-VPS2

Quantity	Component Description
1	Dressing: Tegaderm <sup>®4</sup> 15.5 cm x 10 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: BIOPATCH <sup>®5</sup>
1	Dressing: STATLOCK <sup>®6</sup> PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Surgical Apparel: Mask with Eye Shield
1	Tape: Steri-Strip <sup>®4</sup>
1	HemoHopper <sup>®</sup> Fluid Receptacle

- <sup>1</sup> Licensed under US Patent No. 7,329,412.
- <sup>2</sup> A trademark of Becton, Dickinson and Company.
- <sup>3</sup> A registered trademark of CareFusion or one of its subsidiaries.
- <sup>4</sup> A registered trademark of 3M Company.
- <sup>5</sup> A registered trademark of Johnson & Johnson Corporation.
- <sup>6</sup> A registered trademark of C. R. Bard, Inc.

## Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.