

Representative photo(s): actual product specifications may vary from that shown in the photo

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® VPS TipTracker™ Stylet

PI PICC KIT: 2L 5.5FR 55CM AGBA TIPTRACK

SKU / Article #: CDC-45552-TTS

Highlights

- Antimicrobial: Arrowg+ard Blue Advance*
- Antithrombogenic: Arrowg+ard Blue Advance*
- Sharps Safety
- Pre-loaded with VPS Stylet

Considerations

- *Chlorhexidine
- Not made with natural rubber latex.
- Includes Lidocaine

Sales Quantity/Case: 3

Pressure Injectable Arrowg+ard Blue Advance® Two-Lumen PICC pre-loaded with Arrow® VPS TipTracker™ Stylet

SKU / Article #: CDC-45552-TTS

Quantity	Component Description
1	Two Lumen TaperFree® Catheter with Arrowg+ard Blue Advance® Antimicrobial/Antithrombogenic Protection ¹ : 5.5 Fr. (1.88 mm OD) x 55 cm, Pressure Injectable, Blue FlexTip®, Arrow® TipTracker™ Stylet and T-Port Connector
1	GlideThru™ Peel-Away Sheath: 5.5 Fr. x 2-3/4" (7 cm) Radiopaque over 5.5 Fr. Dilator
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
1	Injection Needle: SafetyGlide™ ² 25 Ga. x 5/8" (1.60 cm)
1	Safety Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW

Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC pre-loaded with Arrow VPS
TipTracker Stylet

SKU / Article #: CDC-45552-TTS


Quantity	Component Description
1	Syringe: 10 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
4	10 mL Luer-Lock Pre-Filled Saline Syringe
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®3} One-Step Solution
1	3 mL Applicator 2% CHG and 70% IPA Chloraprep ^{®3} One-Step Solution with Hi-Lite Orange [™] Tint
2	Dust Cap: Non-Vented
1	SecondSite [™] Adjustable Hub: Fastener
1	SecondSite [™] Adjustable Hub: Catheter Clamp
1	SharpsAway [®] II Locking Disposal Cup
1	SharpsAway [®] Disposal Cup
1	Catheter Trimmer
1	Towel: 24" x 36" (61 cm x 91 cm)
1	Maximal Barrier Drape [™] with two 13 cm (5") x 15 cm (6") fenestrations
1	Towel
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Patient ID Card
1	Chart Sticker
1	Patient Information Booklet
1	Medication Label: 1% Lidocaine
1	Checklist/CLIP Sheet
2	Paper Tape Measure
1	Tourniquet
1	Foam Electrodes (3 per pack)
10	Gauze Pad: 2" x 2" (5 cm x 5 cm)
10	Gauze Pad: 4" x 4" (10 cm x 10 cm)

Pressure Injectable Arrowgard Blue Advance Two-Lumen PICC pre-loaded with Arrow VPS TipTracker Stylet
 SKU / Article #: CDC-45552-TTS

Quantity	Component Description
1	Surgical Apparel: Impervious Gown
1	Dressing: I.V. Advanced 8.5 cm x 11.5 cm
1	Dressing: BIOPATCH ^{®4}
1	Dressing: STATLOCK ^{®5} PICC Plus Catheter Stabilization Device
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask
1	Surgical Apparel: Mask with Eye Shield
1	Equipment Cover: 13 cm x 244 cm (5" x 96")
1	Tape: Steri-Strip ^{®6}
1	HemoHopper [®] Fluid Receptacle
1	T-Piece Cover

- ¹ Licensed under US Patent No. 7,329,412.
- ² A trademark of Becton, Dickinson and Company.
- ³ A registered trademark of CareFusion or one of its subsidiaries.
- ⁴ A registered trademark of Johnson & Johnson Corporation.
- ⁵ A registered trademark of C. R. Bard, Inc.
- ⁶ A registered trademark of 3M Company.

Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

