



Representative photo(s): actual product specifications may vary from that shown in the photo

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

PSI SET: 8 FR

SKU / Article #: CI-09803

Highlights

CE Mark

Considerations


Caution: This product contains dry natural rubber.

Sales Quantity/Case: 10

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters  
SKU / Article #: CI-09803

Quantity	Component Description
1	Arrow-Flex® Sheath: 8 Fr. (2.5 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Obturator: 8 Fr.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.