

Two-Lumen CVC

CVC KIT: 2-LUMEN 7 FR X 16 CM

SKU / Article #: CK-03000-B

Highlights

CE Mark

Considerations




Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC	
SKU / Article #: CK-03000-B	
Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 3 mL Luer-Lock Syringe
1	Syringe: 3 mL Luer-Lock
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
2	Dust Cap: Non-Vented
1	SharpsAway® Disposal Cup
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.