

Arterial Catheterization Set

ARTERIAL CATH SET: 18 GA X 5IN

SKU / Article #: CK-04018

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25


Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set

SKU / Article #: CK-04018

Quantity	Component Description
1	Catheter: 18 Ga. x 12 cm Radiopaque Polyurethane with Integral Suture Wing
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.