



One-Lumen CVC

CVC KIT: 14 GA X 6"

SKU / Article #: CK-04711

Contraindications & Considerations

Not made with natural rubber latex.

Includes Lidocaine

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo



One-Lumen CVC
SKU / Article #: CK-04711

Quantity	Component Description
1	One-Lumen Catheter: 14 Ga. (2.2 mm OD) x 16 cm, Side Holes
1	Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
2	Syringe: 5 cc Luer-Slip
1	Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	3 mL Applicator 2% CHG and 70% IPA ChloraPrep ^{®1} One-Step Solution with Hi-Lite Orange [™] Tint
1	Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
1	Scalpel: #11
1	Dressing: OpSite ^{®2} 7 cm x 6 cm
2	Gauze Pad: 2" x 2"
4	Gauze Pad: 4" x 4"
1	Suture: 3-0 Braided Silk with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

² A registered trademark of T.J. Smith & Nephew, Limited.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

© 2023 Teleflex Incorporated. All Rights Reserved.