

Two-Lumen Hemodialysis Curved Catheter for High Volume Infusions

HEMODIALYSIS SET: 2-LUMEN 14 FR X 15 CM

SKU / Article #: CS-12142-CF

Highlights

CE Mark

Considerations




Not made with natural rubber latex.

Sales Quantity/Case: 5

Two-Lumen Hemodialysis Curved Catheter for High Volume Infusions	
SKU / Article #: CS-12142-CF	
Quantity	Component Description
1	Two-Lumen Catheter, Curved: 14 Fr. (4.72 mm OD) x 15 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
1	Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 16 Fr. (5.3 mm) x 11.4 cm
1	Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
2	Dust Cap: Non-Vented
1	Suture Wing
1	Drape: 24" x 36" (61 cm x 91 cm)
1	Scalpel: #11

Representative photo(s): actual product specifications may vary from that shown in the photo

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.