

Representative photo(s): actual product specifications may vary from that shown in the photo

Pediatric Two-Lumen CVC

CVC SET: 2-LUMEN 4 FR X 13 CM

SKU / Article #: CS-14402

Highlights

CE Mark

Considerations

Not made with natural rubber latex.


Sales Quantity/Case: 5

Pediatric Two-Lumen CVC

SKU / Article #: CS-14402

| Quantity | Component Description |
|----------|---|
| 1 | Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm |
| 1 | Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer |
| 1 | Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle |
| 1 | Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW |
| 1 | Syringe: 5 mL Luer-Slip |
| 1 | Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm |
| 1 | Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm |
| 2 | Dust Cap: Non-Vented |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |
| 1 | SecondSite™ Adjustable Hub: Fastener |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.