

Representative photo(s): actual product specifications may vary from that shown in the photo

Pediatric Two-Lumen CVC

CVC SET: 2L 5 FR X 13 CM

SKU / Article #: CS-14502

Highlights

CE Mark

Considerations

Not made with natural rubber latex.





Sales Quantity/Case: 5

Pediatric Two-Lumen CVC

SKU / Article #: CS-14502

| Quantity | Component Description |
|----------|---|
| 1 | Two-Lumen Catheter: 5 Fr. (1.8 mm OD) x 13 cm |
| 1 | Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) |
| 1 | Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW |
| 1 | Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip Syringe |
| 1 | Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm |
| 1 | Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm |
| 2 | Dust Cap: Non-Vented |
| 1 | SecondSite™ Adjustable Hub: Fastener |
| 1 | SecondSite™ Adjustable Hub: Catheter Clamp |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.